



Uploading Your Omnipod 5 Pump using Glooko

Data is uploaded from the Omnipod 5 Controller wirelessly to Glooko.

Option 1: Using the Glooko App to set up your account

1. Download and Set Up the App:

- Search for the Glooko App in the App Store, Glooko mobile app in the google play store, and download it.
- Open the app and create an account

2. Add Care team

- From the home screen, tap Add Care Teams
- Next screen, tap Add care teams at the bottom of the screen
- Enter proconnect code cbd for Charles H. Best Diabetes Centre
- Tap Continue
- Tap Connect

3. Sync your diabetes device

- From the home screen, tap sync your diabetes device
- Tap continue, follow prompts on the screen to connect your sensor and omnipod 5
- You will be taken to the omnipod 5 website to login with your Omnipod ID and connect

Option 2: Using a Computer to set up your account

There is no option to physically plug in your controller to upload the data to Glooko.

1. Set Up Your Glooko Account:

- Visit Glooko: <https://join.glooko.com/>
- Use proconnect code cbd for Charles H. Best Diabetes Centre

2. Link Glooko to your omnipod ID

- Log into your online omnipod account at www.omnipod.com, using your Omnipod ID
- Go to the My Account tab, select Account Linking from the dropdown menu

- Select log in with Glooko
- Proceed to Glooko and complete the process to link your account

When complete, your status will update to Linked and data will automatically upload into your personal Glooko account.