



Sharing CamAPS Fx data via mylife Cloud

Use mylife Cloud to automatically upload and share data from the CamAPS Fx app.

1. Create mylife cloud account. Go to <https://mylife-software.net/>, click register to create an account. Use the same email address used for the CamAPS Fx app.
2. Open the CamAPS Fx app
3. Open the main menu (3 lines top left), go to account and verify that consent is given to share diabetes data.
4. Go to the main menu, select “share”
5. Under mylife cloud (first section) and user 1, Tap on the blue “none” and log in with your mylife cloud username and password. Press link.
6. Confirm that data synchronization is active by visiting the mylife cloud profile. Note: The initial upload may take up to 48 hours.
7. You will receive an invited from the Best Centre, confirm that you wish to share therapy data. Go to personal data in the mylife cloud account and check the Best Centre to give permission.