



Uploading Your mylife Loop YpsoPump Data


Data is automatically uploaded from the CamAPS FX app wirelessly to glooko.

1. a. Create a personal Glooko account

- Download and Set Up the App: Search for the Glooko mobile app in the google play store, and download it.
- Open the app and create an account

b. Connect Glooko to your clinic

On the homescreen of the Glooko app, tap “Set up your clinic connection”. (You may need to scroll to the right to find this option.)



- Enter proconnect code **cbd** for Charles H. Best Diabetes Centre
- Tap Continue, Connect, Done

2. Connect CamAPS FX app to Glooko

- Open the **Main Menu** ☰
- Click on **Share**
- Look for the **Glooko** section and tap **None** under **User 1**
- Enter your Glooko username and password
- Press **Link**

Glooko support can be found at:

- Call: 1-800-206-6601
- Email: support@glooko.com
- Text: 650-720-5310
- <https://ca.join.glooko.com>