

Uploading Data from Your Freestyle Libre Flash Glucose Monitor

***To save time, PLEASE UPLOAD YOUR DATA BEFORE YOUR APPOINTMENT ***

Email bloodsugars@charleshbest.com or your Diabetes Educator once you have uploaded if you would like the data to be reviewed.

Go to LibreView: <https://www2.libreview.com/>

Reader users:

- Email your educator or bloodsugars@charleshbest.com and ask for invitation to LibreView to share your data
- You will receive an email inviting you to LibreView. Follow the link.
- Select “Sign up” to create an account.
- Accept the privacy and sharing agreements. A pop up will appear asking you to accept sharing invitation from Charles. H. Best Centre, click on “Accept”. You only need to do this once.
- Under “Upload a Device”. Click on “Download the LibreView Device Driver”. Find this file in your downloads folder, open it to install the software.
- Return to LibreView. Plug in your reader with the yellow cord to a USB port on your computer. Click on "upload device".

LibreLink App users:

- Email your educator or bloodsugars@charleshbest.com and ask for invitation to LibreView to share your data.
- You will receive an email inviting you to LibreView. Follow the link and log in using the email and password you used to set up the Libre app on your phone.
- Accept the privacy and sharing agreements. A pop up will appear asking you to accept sharing invitation from Charles. H. Best Centre, click on “Accept”.
- You only need to do this once. After this is completed, to have your data reviewed email bloodsugars@charleshbest.com and ask to look at your Libre data and mention that you are connected to LibreView.

Option 2: Using Practice ID

- Open the Settings Menu: **Click the three blue lines to open the Settings Menu.**
- Click Account Settings

- Click My Practices: **Click on the My Practices section.**
- Link to a practice: Enter the Practice ID in the field provided and click the Add button. The Best Centre clinic ID # is **08878529**. You should now see The Best Centre Clinic name appear in your Linked Practices list. This means you have successfully linked accounts and the Best Centre now has access to all your data.
- *Note: You can link your account to multiple Practices. You can also remove a Practice at any time by clicking the red **Remove** link.*