



## Uploading Your mylife Loop (YpsoPump) Data

Data is automatically uploaded from the CamAPS FX app wirelessly to glooko.

### 1. Create a personal Glooko account

- Download and Set Up the App: Search for the Glooko mobile app in the google play store, and download it.
- Open the app and create an account

### 2. Add Care team

- From the home screen, tap Add Care Teams
- Next screen, tap Add care teams at the bottom of the screen
- Enter proconnect code cbd for Charles H. Best Diabetes Centre
- Tap Continue
- Tap Connect, tap Done

### 3. Connect CamAPS FX app to Glooko

- From the home screen, select drop down menu (top left)
- Tap Share
- Under the Glooko heading, select user 1. Enter Glooko username and password.
- Tap link

### Glooko support can be found at:

- Call: 1-800-206-6601
- Text: 650-720-5310
- Email: [support@glooko.com](mailto:support@glooko.com)

<https://ca.join.glooko.com>