



Uploading Your Medtronic 780g Pump and Guardian Sensor

Option 1: Using the MiniMed Mobile App

1. Download and Set Up the App:

- Search for the MiniMed Mobile App in the Google Play Store or Apple App Store and download it.
- Open the app, select your country. Select “Create account” at the bottom of the screen
- Log in using your Carelink account credentials.
- Follow the app instructions to pair your pump to the app using Bluetooth.

2. Upload Pump Data via the App:

- Open the app and log in.
- Tap the menu (three lines) in the top left corner and select ‘Sync to Carelink.’
- Click ‘Upload Now’ to upload your pump data.

Option 2: Using a Computer and USB Adapter

1. Create a Carelink Account:

- Go to Carelink Minimed EU: <https://carelink.minimed.eu>
- Follow the steps to create an account.

2. Install Uploader Software:

- Log in to your Carelink account and click on ‘How to Install Uploader’ under the ‘Upload Now’ button.
- Download and install the Carelink Uploader. Restart your computer if required.

3. Upload Pump Data:

- Connect the Blue Adapter USB to your computer.
- Log in to Carelink, click ‘Upload Now,’ and follow the instructions to pair your pump.
- Clear any active alarms on the pump and ensure the battery is charged.
- Complete the upload process and click ‘Exit.’