

Please complete the following steps to prepare for your upcoming **MINIMED 780G PUMP UPGRADE**.

1. Open your box(s), take everything out and look for your MiniMed 780G Training guides. Please read these carefully so that you are prepared.
2. Go to <https://www.medtronic.com/ca-en/diabetes/home/support/product-support/minimed-780g-training.html> to access the online learning modules. Please complete the MiniMed 780G training modules prior to your training. Use your pump to follow along with the online tutorial.
3. Place batteries in your pump and familiarize yourself with the various menus and buttons. Enter the time, date and your current settings. These are the settings you need from the pump you are using now:
Insulin: carbohydrate ratios (ICR)
Insulin sensitivity or correction factor (ISF)
Basal programs
Basal rates
BG Targets
Active Insulin time (Insulin action)
Low reservoir
Max basal
Max bolus
Auto off (Auto suspend)
4. Make sure you have ordered, received and charge your Guardian Link Transmitter and ordered and received your sensors. If you require sensors, please order them by calling 1-800-284-4416. Many leave this until the last minute and find that they do not have what they need for training.
5. Go to <http://www.carelink.minimed.eu> and set up your Medtronic CareLink personal account. This will allow you to access the web-based data management tool and share your pump data with your healthcare professional.
6. Download and set up the MiniMed Mobile App onto your smartphone.

You will need the following supplies for your appointment:

- Pump box and pump accessory box with contents
- Infusion sets (please bring at least 2 with you)
- Medtronic pump reservoir (bring the entire box to ensure nothing is missing)
- Blood glucose meter and test strips (even if using a sensor)
- Blood ketone meter and test strips
- Insulin

Please contact Medtronic Technical Support at 1-800-284-4416 if you have any questions as you complete the preparation for training.